0377



Aims

Investigate aldo-keto reductase (AKR) activity in lung cancer cells as an initial step towards developing an **EVOC®** Probe-based breath test for cancer screening and early detection by:

- Evaluating the expression of AKRs associated with aldehyde metabolization in lung cancer and non-neoplastic cancer samples.
- Detecting AKR-associated alcohols and aldehydes from in vitro headspace using gas chromatography-mass spectrometry (GC-MS) workflow as used from clinical breath samples.
- Assessing the impact of AKR inhibition on VOC levels in headspace samples.

1. Background and Objectives

Cancer metabolism represents a promising and largely untapped focus for diagnostic testing. What's more, while the genomics of cancer can vary extensively, these changes converge onto key metabolic pathways bringing about survival benefits that enable cancer cells to survive in the harsh tumour microenvironment. The use of agents such as UDP-glucose in PET scans serves to illustrate how we can target metabolic pathways with exogenous compounds to provide high sensitivity cancer detection.

Our hypothesis is that a similar approach can be used in combination with non-invasive breath sampling to provide reliable early detection of cancer in a form that is well suited to the type of screening programme that could dramatically improve cancer survival.

Lung cancer is an area where there is significant unmet need to improve early detection through screening of at risk populations. Incidence and mortality in these cancers are high

and the majority of cases are diagnosed in the later stages. In addition, the location of lung cancers makes it particularly well suited for detection via breath sampling. As such, we have chosen to make this the initial focus of our work.

Rapid growth, poor blood flow and persistent genetic errors in cancer cells contribute to a high level of oxidative stress characterized by an increase in reactive oxygen species (ROS). In turn, ROS promote destructive processes such as lipid peroxidation which produce aldehydes. Human lung cancers increase expression of aldo-ketoreductase (AKRs) enzymes to help process these excess aldehydes and reduce them to alcohols.

We are seeking to investigate whether AKRs can be targeted using an exogenous volatile organic compound probe (EVOC® Probe), the metabolism of which can be monitored on breath using Breath Biopsy[®].



2. Method



IC50 (µM) Inhibitor AKR1B1 AKR1B10 AKR1B15 0.01 0.006 Tolrestat >100 0.034 ± 0.005 JF0064 0.3 50 µL sample 10 mL headspace Centri Automated Sampling TD-GC-MS Analysis Manual HiSorb PDMS HiSorb BenchToF MS VOC capture Pre-incubat Bench TOP2 30 min 37°C 10 min 37°C agitate 300rpm agitate 300rpm

Figure 2: Method overview of biological and chemical in vitro methods used to investigate the effect of AKR activity modulation on released headspace VOCs. For inhibition of the catalytic activity of AKR1B1/B10/B15, two lung cancer cell lines were treated for 24h with two AKR inhibitors (Tolrestat: IC50 and/or IC75 based on AKR activity assay; JF0064: 10µM). For permanently abrogated AKR1B10 or AKR1B15 expression, three clones of AKR1B10 knockout and three clones of AKR1B15 knockout in the background of A549 cells were developed using CRISPR-Cas9 system. A Mock-Cas9 was also used as control (wild-type cells). Cells were treated with 10μ M of aldehydes (30 µM for Alkanal 1), and aliquots of cell culture media were collected at 3 timepoints for headspace analysis. Samples were stored at -80°C before analysis. An additional plate was run in parallel in each experiment for the control of evaporation and/or cross-contamination. Aldehydes were added to this plate following the plate layout of each experiment set. Positive controls: compound spiked in cell culture media in a closed vial. Blanks: vehicle DMSO.

 Table 1: Inhibitory effect of different compounds on enzymatic
activity in non-cell based assay. D,L-glyceraldehyde as a substrate. JF0064: 2,2',3,3',5,5',6,6'-octafluoro-4,4'- biphenyldiol. Adapted from Giménez-Dejoz et al. (2015; DOI:10.1371/journal.pone.0134506).

Figure 3: Automated and manual headspace VOC sampling methods coupled with a mass spectrometry system. 50 μ L of aqueous samples containing the volatiles are used for headspace analysis. Quality control mix (Deuterated standards and a panel of aldehydes and alcohols) are injected onto sorbent tubes and run throughout the sequence. The mass spectrometry method is selected based on the higher sensitivity and specificity for each pair of aldehydes and alcohol products.TD-GC-MS: Thermal desorption-gas chromatography-mass spectrometry. TOF-MS: Time-of-flight mass spectrometry.

Using exogenous volatile organic compound (EVOC) probes to target tumour-associated aldo-keto reductase activity: a potential tool to detect lung cancer

3. Results and Discussion



Figure 4: Expression of AKRs in samples from lung cancer patients. (A) AKR1B1, AKR1B10 and AKR1B15 expression in lung cancer samples and non-neoplastic samples (cancer free surgical margin) from patient with lung cancer from TCGA database; most of cancer samples discriminate from non-neoplastic samples based on these three AKR gene expression. (B) Representative immunostaining of a lung cancer specimens showing high expression of AKR1B1 and/or AKR1B10 in lung cancer cells than in the surrounding non-neoplastic lung cells. TCGA: The Cancer Genome Atlas (TCGA) (https://www.cancer.gov/ccg/research/genome-sequencing/tcga). TCGA-LUAD: Lung adenocarcinomal dataset. TCGA-LUSC: Lung Squamous Cell Carcinoma.



Figure 5: Example of calibration curves for AKR alkanal substrate and the alkanol bioproduct. High sensitivity and broad linearity can be observed.

A549 CELLS

H460 Cells Cells vs Cell Culture Media Only





Figure 7: AKR inhibitor effects in lung cancer cells. A commercial colorimetric AKR activity assay validated to AKR1B10, AKR1C1 and AKR1C3 activity analysis was used to evaluate the effect of Tolrestat and JF0064 in two lung cancer cell lines. A) Dose-dependent effect of Tolrestat in AKR activity in A549 and H460 after 24h of treatment. B) JF0064 does not affect overall AKR activity after 24h of treatment in A549 and H460 cell lines.

Figure 1. The EVOC Probe approach. An administered EVOC Probe is introduced and monitored non-invasively on breath along with its metabolic bioproducts. Aldo-keto reductases (AKRs) are a potential EVOC Probe target that are upregulated in some cancers as an adaptation to oxidative stress and lipid peroxidation.



Mariana Ferreira Leal¹, Alexandra. Martin², Ben Taylor², Rory Stallard², Connor Clarke², Christiaan F Labuschagne¹, Rob Smith², Matthew Hart³, Max Allsworth³, Billy Boyle³ ¹Owlstone Medical, Clinical and Translational Science, Cambridge, United Kingdom, ²Owlstone Medical, Analytical Science, Cambridge, United Kingdom, ³Owlstone Medical, Leadership team, Cambridge, United Kingdom. *email: breathbiopsy@owlstone.co.uk

ostrate	LOD (µM)	LOQ (µM)	Bioproduct	LOD (µM)	LOQ (µM)
kanal 1	0.840	2.546	Alkanol 1	0.307	0.931
anal 2	0.14	0.43	Alkanol 2	0.04	0.43
kenal 1	0.03	0.08	Alkenol 1	0.15	0.45
enal 2	0.05	0.14	Alkenol 2	0.15	0.46
c Aldehyde 1	0.07	0.22	Aromatic Alcohol 1	0.16	0.48

Table 2: Sensitivity of analytical method to detect aldehyde substrates and alcohol bioproducts using high-capacity PDMS sorbent probes.



Figure 6: Monitoring AKR activity in lung cancer cells using aldehydes as EVOC probes. Aldehydes were added to cell culture media with and without lung cancer cells (evaporation controls) in 24-well plates. Headspace from A549 and H460 samples showed lower levels of aldehydes and higher levels of alcohol products than those observed in evaporation controls confirming that AKRs are active in the studied lung cancer cells and can be monitored using the headspace system. VOC: volatile organic compounds (aldehydes and alcohols)







4. Conclusions

A549 mock-Cas9 cells.

• AKR enzymes are potential targets for EVOC probes to detect lung cancer.

• We have demonstrated the capacity to detect VOCs from in vitro headspace. Using in vitro study of human lung cancer cells, we have demonstrated the potential to monitor the metabolic conversion of administered EVOC Probe aldehydes into alcohols by AKRs, and have used AKR inhibition/silencing to show that this relationship is specific and sensitive to manipulation.

• Pre-clinical samples were analysed using the same GC-MS workflow used for the detection of VOCs in clinical breath samples (Breath Biopsy workflow). We are now seeking to investigate the same relationships through in vivo and ex vivo sampling to further establish the potential to use EVOC Probes targeting AKR metabolism as a tool for early detection of lung cancer.

