

BREATH BIOPSY

Dynamic Limonene Breath Testing Maximizes Classification Performance for Subjects with Cirrhosis

^{1,8}Ferrandino G. ¹Ricciardi F. ¹Murgia, A. ¹Banda I. ¹Manhota M. ¹Ahmed Y. ¹Sweeney K. ¹Nicholson-Scott L. ¹McConville L. ¹Gandelman O. ¹Allsworth, M. ¹Boyle B. ^{1,2}Smolinska A. ³Ginesta A. ³Contreras J. ⁴Asenjo-Lobos C. ⁵Barrientos V. ⁵Clavo N. ⁶Novoa A. ⁴Riviotta A. ⁷Jerez M. ^{3,8}Méndez L.

¹ Owlstone Medical, 183 Cambridge Science Park, Milton Road, Cambridge CB4 0GJ, UK, Email: breathbiopsy@owlstone.co.uk
² Department of Pharmacology and Toxicology, School for Nutrition and Translational Research in Metabolism (NUTRIM), Maastricht University Medical Center, the Netherlands.
³ Gastroenterología, Clínica Alemana, Facultad de Medicina Clínica Alemana Universidad de Desarrollo, Santiago, Chile
⁴ Centro de Estudios Clínicos, Instituto de Ciencias e Innovación en Medicina (ICIM), Facultad de Medicina Clínica Alemana Universidad del Desarrollo, Santiago, Chile
⁵ Unidad de endoscopia, Hospital Padre Hurtado, Santiago, Chile.
⁶ Laboratorio de Fisiología Digestiva, Clínica Alemana, Santiago, Chile.
⁷ Nursing student, Universidad de Las Américas, Santiago, Chile.
⁸ Corresponding author

- Aims
- Assess diagnostic performance of a dynamic limonene breath test for cirrhosis detection.
 - Investigate limonene exhalation kinetics after oral administration.
 - Explore the correlation of limonene bioavailability with cirrhosis severity.

1. Background and Objectives




Figure 1: Progression of chronic liver disease to cirrhosis is often asymptomatic. Patients that are diagnosed at earlier stages benefit from causal treatments and can achieve a life expectancy similar to that of the healthy population. On the contrary, patients diagnosed with episodes of decompensation (liver failure) have a poor prognosis. Methods to detect cirrhosis in primary care suffer from low performance and most tests performed are relegated to secondary or tertiary care. This short diagnostic pathway results in reports of up to 75% of patients receiving a diagnosis of cirrhosis only upon decompensation¹.




Figure 2: Altered levels of limonene alongside certain additional compounds in exhaled breath can identify subjects with cirrhosis compared to controls, demonstrating high diagnostic performance. However, limonene and the other compounds originate from the diet therefore the amount of these compounds in the breath depends on the subject's diet^{2,3}. Therefore, diet is a confounding factor that needs to be controlled to maximize the diagnostic performance of compounds in exhaled breath for cirrhosis.

References

1. Labenz C, Arslanow A, Nguyen-Tat M, Nagel M, Worms MA, Reichert MC, et al. Structured Early detection of Asymptomatic Liver Cirrhosis: Results of the population-based liver screening program SEAL. *J Hepatol* 2022;77(3):695-701. doi: 10.1016/j.jhep.2022.04.009. PMID: 35472313

2. Ferrandino G, Orf I, Smith R, Calcagno M, Thind AK, Debram-Beecham I, et al. Breath Biopsy Assessment of Liver Disease Using an Exogenous Volatile Organic Compound-Toward Improved Detection of Liver Impairment. *Clin Transl Gastroenterol* 2020;11(9):e00239. doi: 10.14309/ctg.0000000000000239. PMID: 33094960

3. Ferrandino G, De Palo G, Murgia A, Birch O, Tawfike A, Smith R, et al. Breath Biopsy((R)) to Identify Exhaled Volatile Organic Compounds Biomarkers for Liver Cirrhosis Detection. *J Clin Transl Hepatol* 2023;11(3):638-648. doi: 10.14218/JCTH.2022.00309. PMID: 36969895

4. Fernandez Del Rio R, O'Hara ME, Holt A, Pemberton P, Shah T, Whitehouse T, et al. Volatile Biomarkers in Breath Associated With Liver Cirrhosis - Comparisons of Pre- and Post-liver Transplant Breath Samples. *EBioMedicine* 2015;2(9):1243-1250. doi: 10.1016/j.ebiom.2015.07.027. PMID: 26501124

2. Methods




Figure 3: This cross-sectional case-control study was approved by “Comité Ético Científico de la Facultad de Medicina Clínica Alemana Universidad del Desarrollo” and all participants provided written informed consent. All of the procedures were conducted in compliance with applicable guidelines for the ethical conduct of the study with their origins in the Declaration of Helsinki. Presence or absence of cirrhosis was established using ultrasound according to EASL and AASLD guidelines. Subjects older than 18 years old and weighing more than 60 kg were randomly enrolled between January 2022 and December 2022 and instructed before the experiment to fast for ≥10 h, to not consume alcohol the previous day, and to not brush their teeth or use mouthwash in the previous 2 h. Breath Biopsy® samples were collected before and after ingestion of 100 mg limonene. Samples were shipped to Owlstone Medical (UK) for analysis using gas-chromatography mass spectrometry (GC-MS). Absolute limonene quantification was obtained against a standard curve.

3. Results




Table 2. Limonene exhalation kinetic parameters			
Parameter	Control	Cirrhosis	p-value
C _{max} (ng)	595	2077	< 0.001
Log ₁₀ C ₀ median [IQR]	6.9 [361 - 903]	8.4 [1051 - 4260]	< 0.001
T _{max} n (%)	18 (62.1%) 40 min (37.9%)	13 (44.8%) 12 (41.4%) 2 (6.9%) 0	-
AUC (0-90 min) ng*min/400ml median [IQR]	27107 [17605 - 34946]	121437 [57921 - 202733]	<0.001
Slope	-0.027 [-0.031; -0.023]	-0.025 [-0.027 -0.019]	0.072

Figure 4: After overnight fasting, subjects with cirrhosis showed significantly elevated levels of breath limonene as expected⁴ (A). Limonene administration induced a spike on breath with significantly higher levels in subjects with cirrhosis at all of the tested timepoints (B). A higher limonene bioavailability was observed in the cirrhosis group (C). A semi logarithmic presentation of exhaled limonene profile showed a first order decay with > 90% of the subjects showing an R² > 0.8 (D). The slope of the control and cirrhosis group showed no significant difference (E) while the cirrhosis group showed higher c₀ (intercept with y axes) (F).






Figure 5: Levels of limonene on breath at baseline after an overnight fast showed high diagnostic/ classification performance for subjects with cirrhosis as reported previously⁴ (A). All of the timepoints post administration showed improved classification performance compared to baseline, with the best performing time point being at 60 minutes (B). Importantly, limonene administration boosted sensitivity, the main parameter affected at baseline (C).

Table 3. Summary of limonene diagnostic performances at different timepoints				
Timepoint (min)	AUROC	Sensitivity/ Specificity (%)	+/- Predictive Values (%)	+/- Likelihood Values
0	0.83 ± 0.12	0.66 ± 0.09/ 0.83 ± 0.07	79.17/70.59	3.8/0.42
20	0.92 ± 0.07	0.82 ± 0.095/ 0.79 ± 0.08	79.88/81.62	3.97/0.23
40	0.94 ± 0.06	0.79 ± 0.08/ 0.9 ± 0.06	88.37/80.71	7.6/0.24
60	0.91 ± 0.07	0.83 ± 0.07/ 0.9 ± 0.06	88.89/83.87	8.0/0.19
90	0.91 ± 0.07	0.76 ± 0.08/ 0.86 ± 0.06	84.62/78.12	5.5/0.28
120	0.93 ± 0.06	0.83 ± 0.07/ 0.86 ± 0.06	85.71/83.33	6.0/0.2



4. Conclusions and Next Steps

- Dynamic breath analysis enhances classification performance for cirrhosis opens up the opportunity to implement a reliable test in primary care.
- Alterations of limonene breath kinetics in cirrhosis resemble the changes induced by drugs with high hepatic extraction.
- Correlation of limonene with disease severity suggests applications for non-invasive monitoring of therapeutic interventions.
- Identify additional compounds with low hepatic extraction to assess their accuracy in reporting on metabolic alterations.
- Explore applications for disease progression and regression.
- Investigate similar approaches in earlier stages of liver disease such as non-alcoholic steatohepatitis (NASH).