BREATH BIOPSY

Aims

- Assess diagnostic performance of a dynamic limonene breath test for cirrhosis detection.
- Investigate limonene exhalation kinetics after oral administration.
- Explore the correlation of limonene bioavailability with cirrhosis severity.

1. Background and Objectives



receiving a diagnosis of cirrhosis only upon decompensation¹.



Figure 2: Altered levels of limonene alongside certain additional compounds in exhaled breath can identify subjects with cirrhosis compared to controls, demonstrating high diagnostic performance. However, limonene and the other compounds originate from the diet, therefore the amount of these compounds in the breath depends on the subject's diet^{2,3}. Therefore, diet is a confounding factor that needs to be controlled to maximize the diagnostic performance of compounds in exhaled breath for cirrhosis.

References

- 1. Labenz C, Arslanow A, Nguyen-Tat M, Nagel M, Worns MA, Reichert MC, et al. Structured Early detection of Asymptomatic Liver Cirrhosis: Results of the population-based liver screening program SEAL. J Hepatol 2022;77(3):695-701. doi: 10.1016/j.jhep.2022.04.009. PMID: 35472313
- 2. Ferrandino G, Orf I, Smith R, Calcagno M, Thind AK, Debiram-Beecham I, et al. Breath Biopsy Assessment of Liver Disease Using an Exogenous Volatile Organic Compound-Toward Improved Detection of Liver Impairment. Clin Transl Gastroenterol 2020;11(9):e00239. doi: 10.14309/ctg.0000000000000239. PMID: 33094960
- 3. Ferrandino G, De Palo G, Murgia A, Birch O, Tawfike A, Smith R, et al. Breath Biopsy((R)) to Identify Exhaled Volatile Organic Compounds Biomarkers for Liver Cirrhosis Detection. J Clin Transl Hepatol 2023;11(3):638-648. doi: 10.14218/JCTH.2022.00309. PMID: 36969895
- 4. Fernandez Del Rio R, O'Hara ME, Holt A, Pemberton P, Shah T, Whitehouse T, et al. Volatile Biomarkers in Breath Associated With Liver Cirrhosis Comparisons of Pre- and Post-liver Transplant Breath Samples. EBioMedicine 2015;2(9):1243-1250. doi: 10.1016/j.ebiom.2015.07.027. PMID: 26501124

breathbiopsy@owlstone.co.uk

Dynamic Limonene Breath Testing Maximizes Classification Performance for Subjects with Cirrhosis

¹ Owlstone Medical, 183 Cambridge Science Park, Milton Road, Cambridge CB4 OGJ, UK, Email: breathbiopsy@owlstone.co.uk ² Department of Pharmacology and Toxicology, School for Nutrition and Translational Research in Metabolism (NUTRIM), Maastricht University Medical Center, the Netherlands. ³ Gastroenterología, Clínica Alemana, Facultad de Medicina Clínica Alemana Universidad de Desarrollo, Santiago, Chile ⁴ Centro de Estudios Clínicos, Instituto de Ciencias e Innovación en Medicina (ICIM), Facultad de Medicina Clínica Alemana Universidad del Desarrollo, Santiago, Chile ⁵ Unidad de endoscopia, Hospital Padre Hurtado, Santiago, Chile. ⁶ Laboratorio de Fisiología Digestiva, Clínica Alemana, Santiago, Chile. ⁷ Nursing student, Universidad de Las Américas, Santiago, Chile. ⁸ Corresponding author

Up to 75% of the cases detected in secondary care with decompensation Cirrhosis vs. Healthy (Training Set) Mean ROC (AUROC = 0.99) ±1 std. dev (0.00) Cirrhosis vs. Healthy (Test Set) Mean ROC (AUROC = 0.95) ±1 std. dev (0.04)

2. Methods



Figure 3: This cross-sectional case-control study was approved by "Comité Ético Científico de la Facultad de Medicina Clínica Alemana Universidad del Desarollo" and all participants provided written informed consent. All of the procedures were conducted in compliance with applicable guidelines for the ethical conduct of the study with their origins in the Declaration of Helsinki. Presence or absence of cirrhosis was established using ultrasound according to EASL and AASLD guidelines. Subjects older than 18 years old and weighing more than 60 kg were randomly enrolled between January 2022 and December 2022 and instructed before the experiment to fast for ≥ 10 h, to not consume alcohol the previous day, and to not brush their teeth or use mouthwash in the previous 2 h. Breath Biopsy[®] samples were collected before and after ingestion of 100 mg limonene. Samples were shipped to Owlstone Medical (UK) for analysis using gas-chromatography mass spectrometry (GC-MS). Absolute limonene quantification was obtained against a standard curve.

3. Results

Table 1	Control	Cirrhosis	p-values	А
Number of patients	29	29	-	120 -
Age median [IQR] years	43 [38-57]	59 [54-67]	< 0.001	80 -
Sex n (%)	M - 11 (38%) F - 18 (62%)	M - 9 (31%) F - 20 (69%)	-	0 - 00 -
Height median [IQR] cm	163 [158-170]	160 [156-170]	0.32	40 -
Weight median [IQR] Kg	75 [64-84]	78 [68-85]	0.50	20 -
BMI median [IQR]	26.8 [23.9-31.2]	27.7 [26.0-32.8]	0.23	0 -
Waist circumference median [IQR] (cm)	90 [81.5-104]	106 [97-113]	0.036	C
Child-Pugh class	-	23 (78%), 5 (16%), 2 (6%)	-	2500
MELD median [IQR]	-	10 [7.2-12.8]	-	e 1500
FIB4 median [IQR]	1.4 [0.8-3.1]	23.1 [1.9-4]	p < 0.001	1000 L
APRI median [IQR]	0.2 [0.2-0.3]	0.6 [0.4-0.9]	p < 0.001	500
Platelets median [IQR] *10 ⁹ /L	240 [216-294]	147 [107-209]	p < 0.001	0
Total bilirubin median [IQR] (μmol/L)	8.2 [6.5-12.9]	14.3 [8.2-18.4]	p < 0.001	E
Serum Albumin median [IQR] (g/L)	45 [44-45]	40 [37-44.5]	p < 0.001	
INR median [IQR]	1 [1-1.06]	1.2 [1.0-1.4]	p < 0.001	-0.02
ALT median [IQR] (IU/L)	18 [15-26]	28.5 [18.7-38]	p < 0.001	-0.03
AST median [IQR] (IU/L)	21 [18.5-24]	35 [27.2-45.7]	p < 0.001	Siope Siope
GGT median [IQR] (IU/L)	24 [16-32]	76 [59.2-108.7]	p < 0.001	-0.04
ALP median [IQR] (IU/L)	85 [71.5-101.5]	117 [94.5-153]	p < 0.001	-0.05
Creatinine median [IQR] (mg/dL)	0.78 [0.68-1.93]	0.74 [0.64-0.87]	0.31	
Sodium median	141 [138,7-142.5]	142 [139.5-142.5]	0,65	

Table 2. Limonene exhalation kinetic parameters						
Parameter	Control	Cirrhosis	p-value			
Cmax (ng) median [IQR]	595 [361 - 903]	2077 [1051 - 4260]	< 0.001			
Log10 C _o median [IQR]	6.9 [6.69 - 7.29]	8.4 [7.9-8.9]	< 0.001			
Tmax n (%) 20 min 40 min 60 min 90 min 120 min	18 (62.1%) 11 (37.9%) 0 0 0	13 (44.8%) 12 (41.4%) 2 (6.9%) 1 (3.4%) 1 (3.4%)	_			
AUC (0-90 min) ng*min/400ml median [IQR]	27107 [17605 - 34946]	121437 [57921 - 202733]	<0.001			
Slope	-0.027 [-0.031; -0.023;]	-0.025 [-0.027 -0.019]	0.072			

^{1,8}Ferrandino G. ¹Ricciardi F. ¹Murgia, A. ¹Banda I. ¹Manhota M. ¹Ahmed Y. ¹Sweeney K. ¹Nicholson-Scott L. ¹McConville L. ¹Gandelman O. ¹Allsworth, M. ¹Boyle B. ^{1,2}Smolinska A. ³Ginesta A. ³Contreras J. ⁴Asenjo-Lobos C. ⁵Barrientos V. ⁵Clavo N. ⁵Gabrielli V. ⁶Novoa A. ⁴Riviotta A. ⁷Jerez M. ^{3,5}Méndez L.

at 20, 40, 60, 90, 120 min

with GC-MS



Figure 4: After overnight fasting, subjects with cirrhosis showed significantly elevated levels of breath limonene as expected⁴ (A). Limonene administration induced a spike on breath with significantly higher levels in subjects with cirrhosis at all of the tested timepoints (B). A higher limonene bioavailability was observed in the cirrhosis group (C). A semi logarithmic presentation of exhaled limonene profile showed a first order decay with > 90% of the subjects showing an R2 > 0.8 (D). The slope of the control and cirrhosis group showed no significant difference (E) while the cirrhosis group showed higher cO (intercept with y axes) (F).



Figure 5: Levels of limonene on breath at baseline after an overnight fast showed high diagnostic/ classification performance for subjects with cirrhosis as reported previously⁴ (A). All of the timepoints post administration showed improved classification performance compared to baseline, with the best performing time point being at 60 minutes (B). Importantly, limonene administration boosted sensitivity, the main parameter affected at baseline (C).





Figure 6: Limonene bioavailability measured as area under the curve (AUC) in subjects with cirrhosis correlates with disease severity estimated using MELD, and risk of advanced fibrosis estimated using FIB4 and APRI. Canonical correspondence analysis (CCA) score plot using the first canonical variates of limonene bioavailability, estimated as AUC, and scoring parameters. Each projected data point represents the combined information of limonene AUC and scoring systems of one cirrhotic patient (A). Canonical loadings of the scoring systems representing their correlation between the original variable and its canonical variate and express the contribution of each variable to the overall correlation (B). Correlation matrix visualizing the Spearman correlation between indicated variables (C). Intensity of the red and blue colors indicate strength of respectively positive and negative correlations, while circle size represents the significance test between the variables. A total of six subjects with incomplete metadata were excluded from analysis.



4. Conclusions and Next Steps

- reliable test in primary care.
- interventions.

Table 3. Summary of limonene diagnostic performances at different timepoints							
(min)	AUROC	Sensitivity/ Specificity (%)	+/- Predictive Values (%)	+/- Likelihood Values (%)			
	0.83 ± 0.12	0.66 ± 0.09/ 0.83 ± 0.07	79.17/70.59	3.8/0.42			
	0.92 ± 0.07	0.82 ± 0.095/ 0.79 ± 0.08	79.88/81.62	3.97/0.23			
	0.94 ± 0.06	0.79 ± 0.08/ 0.9 ± 0.06	88.37/80.71	7.6/0.24			
	0.91 ± 0.07	0.83 ± 0.07/ 0.9 ± 0.06	88.89/83.87	8.0/0.19			
	0.91 ± 0.07	0.76 ± 0.08/ 0.86 ± 0.06	84.62/78.12	5.5/0.28			
	0.93 ± 0.06	0.83 ± 0.07/ 0.86 ± 0.06	85.71/83.33	6.0/0.2			

Figure 7: Subjects found to have been allocated in the wrong group after ultrasound confirmation of liver condition. Ultrasound image of a subject who received a diagnosis of cirrhosis by enrolling in this study, and initially allocated as healthy (A). Ultrasound image of a subject diagnosed with cirrhosis three years before the breath tests, who showed recovery after treatment (B). Breath profile of these subjects. Shaded areas represent the 95% confidence interval for the cirrhosis group (blue) and control (grey) for comparison. Each line represents a subject (C).

• Dynamic breath analysis enhances classification performance for cirrhosis opens up the opportunity to implement a

• Alterations of limonene breath kinetics in cirrhosis resemble the changes induced by drugs with high hepatic extraction. • Correlation of limonene with disease severity suggests applications for non-invasive monitoring of therapeutic

 Identify additional compounds with low hepatic extraction to assess their accuracy in reporting on metabolic alterations. Explore applications for disease progression and regression. Investigate similar approaches in earlier stages of liver disease such as non-alcoholic steatohepatitis (NASH).