

Introduction

Exhaled breath contains thousands of Volatile Organic Compounds (VOCs), which are products of metabolic activity and biomarkers for diseases. These can be from endogenous processes in the body and exogenous sources such as environmental exposures and pharmaceuticals.


Reliable longitudinal evaluation of breath VOCs are an essential step enabling studies into:


- Understanding biological variability
- Enabling disease monitoring
- Measuring the response to therapeutic interventions
- Assessing effects of environmental exposures
- Evaluating pharmacokinetics

This study aims to examine the long term stability of the breath biopsy platform as well as its ability to detect changes over time relative to biological variability.

To evaluate this, two experiments were conducted:

- Analytical replicates: Analysis of a target set of VOCs in breath following ingestion of a peppermint capsule and the resulting "washout curve" for target over 8 hours (Figure 1) collected as 4 concurrent samples per timepoint.
- Biological variability: 9 repeats of the same experiment over a 5 weeks period in the same individual, sampled at three time points.





Repeated breath collects Replicate VOC sample tubes collected and analyzed at each timepoint

Methods

Breath Sampling

The study subject was a 32-year old, healthy male who standardised diet and exercise for the duration of the study. Baseline collects were performed at two time points. Thereafter the subject ingested one Obbekjaers 200mg peppermint pill. Samples were collected at 16 time points in the initial study. In the longitudinal study timepoints estimating baseline, peak and washout period were selected at -60, 45, and 180 min. Breath samples were collected using the ReCIVA Breath Sampler and analyzed using the Breath Biopsy platform in the Breath Biopsy Clinical Laboratory (Owlstone Medical Ltd, UK).

Sample Analysis and Data-processing

- Analytical replicates: Analysis of a target set of VOCs in breath following ingestion of a peppermint capsule and the resulting "washout curve" for target over 8 hours (Figure 1) collected as 4 concurrent samples per timepoint.
- Transfer onto column using split injection (Agilent Technologies)
- GC-MS analysis utilizing BENCHTOF HD mass spectrometer
- Peak area extraction using MassHunter Quant based on ion signature

Validation of a methodology for evaluating longitudinal change of VOCs in breath

J. Boschmans (1), R. Smith (1), R. Parris (1), B. Boyle (1), D. Apthorp (1), S. Kitchen (1), M. Allsworth, M. P. C. van der Schee (1)

(1) Owlstone Medical Ltd., 183 Cambridge Science Park, Cambridge, CB4 0GJ, UK. marc.vanderschee@owlstone.co.uk

Breath Collection





Figure 2. The ReCIVA Breath Sampler enables reliable, reproducible collection of breath VOCs and pre-concentration for enhanced sensitivity. Pressure and CO2 sensors in ReCIVA provide real-time monitoring of the patient's breathing, allowing different breath fractions to be sampled in a single collection event if desired.

Result: Washout Curve Over 8 Hours







Figure 3. Washout curve following ingestion of a peppermint capsule. Breath samples from repeated collects from one individual over 8 hours (16 timepoints) were analyzed. Four replicate VOC sample tubes were collected and analyzed at each time point.

owlstonemedical.com

Result: Longitudinal Study Over 5 Weeks



collected and analyzed at each time point.

Conclusions

- in a reliable and reproducible way.

References

van der Schee, M. P. *et al.* Breathomics in lung disease. Chest 147, 224–31 (2015) L. M. Heaney et al. Real-time monitoring of exhaled volatiles using atmospheric pressure chemical ionization on a compact mass spectrometer. Bioanalysis 8 (13), 1325-1336 (2016). J. Beauchamp et al. Real-time breath gas analysis for pharmacokinetics: monitoring exhaled breath by on-line proton-transfer-reaction mass spectrometry after ingestion of eucalyptol-containing capsules. J. Breath Res. 4 (2), 026006 (2010). L. Schmidt. Human metabolism of α -pinene and metabolite kinetics after oral administration. Arch. Toxicol. 91 (2), 677-687 (2017). M. M. L. Steeghs. Collision induced dissociation study of 10 monoterpenes for identification in trace gas measurements using the newly developed proton-transfer reaction ion trap mass spectrometer. Int. J. Mass. Spectrom. 263 (2-3), 204-212 (2007) L. M. Heaney. Exhaled breath analysis in exercise and health. A Doctoral Thesis submitted in partial fulfilment of the requirements for the award of Doctor of Philosophy of Loughborough University (UK)(2016).



Figure 4. Longitudinal study repeating the washout experiment nine times over five weeks in a single individual. Breath samples from pre-ingestion control, and time points 1 and 2 (T1 and T2) were analyzed in each washout experiment. Two replicate VOC sample tubes were

The study demonstrates that the Breath Biopsy platform can be used to study longitudinal changes of exhaled VOCs

The evaluation of a washout curve of peppermint oil components provides a model for the study of absorption, distribution, metabolism and excretion of pharmaceuticals.