

Breath Diagnostics

Exhaled breath contains a vast untapped resource of chemical information represented by hundreds of volatile organic compounds. As products of underlying metabolic activity, they directly reflect the state of cells and tissues and are fast becoming recognised as valuable disease biomarkers

By Dr Marc van der Schee
and Dr Theodore W Wilson
at Owlstone Medical

The metabolic information provided by volatile organic compound (VOC) biomarker analysis offers possibilities beyond those generated by genomic data, which may provide a blueprint for disease predisposition but can lack detail on actual disease activity. Coupled with the non-invasive nature of breath sampling, using this method to perform a 'breath biopsy' is a promising approach for a wide range of diseases.

A growing body of evidence shows the potential of VOCs for directing treatment of disease. Their clinical application has now been accelerated by the development of new technologies for the reliable collection of patient breath samples and advances in VOC analysis technology – both of which have the potential to bring breath biopsy to point-of-care use when combined.

VOC Biomarkers in Breath

Breath is particularly valuable as a source of disease biomarkers for respiratory diseases, because it contains VOCs originating from local airways tissues. Biomarkers in breath are by no means limited to pulmonary conditions – breath also contains systemic VOCs from the circulatory system providing

information about the rest of the body, as well as VOCs from exogenous sources (see Figure 1). Hence, VOC biomarkers in exhaled breath represent a rapid and non-invasive means to gain a detailed insight into the functioning of the whole body.

Using the composition of exhaled breath to make clinical decisions is not a new idea – knowledge of distinctive breath odours related to specific diseases stretches back to Hippocrates. Despite this historical pedigree and the obvious promise held by breath analysis, very few breath-based medical tests have found a place in doctor's toolboxes. This is largely due to a lack of standardised breath sampling technology that has resulted in poor repeatability between studies looking for disease-specific VOC biomarkers (2,3).

Breath analysis research usually employs high-end analytical instrumentation such as gas chromatography mass spectrometry (GC-MS) and selected-ion flow-tube MS. These instruments are the gold standard, but are also expensive, bulky and require expertise to operate correctly – hampering both their deployment in clinical scenarios and the transition of VOCs from research laboratories to deployable clinical tests.

New Technology

Recent advances in breath sampling and analysis technology are now

addressing the field's long-standing issues and herald the arrival of an era in which breath biopsy will be as ubiquitous as blood and urine tests. Arguably, the crucial step of technological standardisation in breath sampling has now been taken with the development of a breath sampler designed by the Breathe Free Consortium, a community of more than 100 breath researchers (see Figure 2) (4).

Pressure and CO₂ sensors provide real time monitoring of patient breathing, allowing specific breath fractions to be selected. The device enables different fractions such as bronchial or end-tidal to be stored independently on sorbent tubes, which can then be analysed immediately or saved for later use. Sorbent tubes ensure reproducible collection of VOCs, unlike other commonly used collection media such as Tedlar® bags (5). The sampler is also designed to incorporate a clean air supply in order to prevent interference from unwanted local VOC contaminants.

Another key advance in breath analysis is the miniaturisation of technologies which have the potential to be deployed at the point-of-care. 'Electronic Nose' (eNose) systems have been developed, using arrays of non-selective chemical sensors to detect VOCs. However, they often suffer from unwanted chemical interactions between VOCs and the sensors, low intra-device repeatability, limited temporal stability and poor chemical selectivity (6).

Keywords

Asthma
Biomarkers
Breath
Volatile organic compounds




Figure 1: Schematic representation of alveolus detailing origins of VOCs.
 1. Local cellular metabolites
 2. Exogeneous compounds, eg microbiota
 3. Compounds originating from elsewhere in the body (1)

Field asymmetric ion mobility spectrometry (FAIMS) is now available in a chip-based miniaturised form that offers improved stability and sensitivity when compared to eNose systems. FAIMS distinguishes ionised VOCs according to differences in the speed that they move through a buffer gas under the influence of an oscillating electric field (7). The FAIMS chip is also programmable in software, making it suitable for a wide range of medical applications (6). Along with the Breathe Free sampler, FAIMS is being incorporated into large clinical trials relating to lung cancer, colorectal cancer and asthma (8-10).

Asthma: A Case Study

As many as 334 million people have asthma and up to 14% of the world's children experience its symptoms (11). Asthma management is focused on achieving control of symptoms to minimise the risk of future exacerbations. However, a lot of patients do not respond sufficiently to treatment, so their symptoms are not kept adequately




Figure 2: The Breathe Free open source breath sampler, designed by over 100 collaborating breath scientists

under control. In the UK, 4.4% of patients fail to respond to standard therapies even at high doses, and account for more than 50% of asthma costs to the NHS.

Asthma symptoms are caused by chronic inflammation in the small airways. Many underlying mechanisms are responsible for this inflammation, with different treatments required for different inflammation types. There is no single diagnostic test for asthma at present, in part because of its heterogeneous underlying pathophysiology.





Figure 3: A typical FAIMS spectrum for exhaled breath (7)

The difference in cell activity related to different inflammatory subtypes is reflected by the metabolites produced by these cells, and consequently by the VOCs in exhaled breath. This is why asthma

Figure 4: Predicting steroid responsiveness in steroid-naïve patients (18)



Source: van der Schee et al, CEA, 2013

has become one of the most studied diseases with regard to breath VOCs, and improvements in technology are accelerating progress in the field.

Improved Diagnostics

One breath-based biomarker, fractional exhaled nitric oxide (F_eNO), is already being used to support asthma diagnosis. It is useful because it reflects a protective biochemical pathway of the lungs, but it does lack specificity as it is a single biomarker affected by many processes other than asthma. A study of VOCs in asthma found that they outperform F_eNO and lung function tests in discriminating asthmatics and healthy controls. Furthermore, combining VOCs with F_eNO gave an even higher diagnostic performance (12).

It has also been shown that VOCs can be used to discriminate with a high degree of accuracy between asthmatic children and those with transient wheezing (13). In both this and two other studies, increases in methylated alkanes were associated with asthma prevalence (14,15). These VOCs are potential markers of oxidative stress that have been linked to airway inflammation.

Data collected so far in one clinical trial indicate that VOC biomarkers collected using the Breathe Free sampler and measured using GC-MS combined with FAIMS can differentiate patients with asthma from patients with a wide range of other pulmonary diseases with a high degree of accuracy (8).

The existing evidence suggests that VOCs are strongly affected by the airway inflammation that characterises obstructive pulmonary diseases like asthma, and that they could be useful tools in non-invasive asthma diagnosis. This method is particularly valuable for diagnosis in children as other known biomarkers – F_eNO and blood immunoglobulin E (IgE) levels, for example – require exhalation at a continuous flow or invasive sampling.

Asthma Phenotyping

While the diagnostic power of breath VOCs could provide a valuable new tool for clinicians, the most impactful application of an asthma breath biopsy is likely to be in the phenotyping of the disease by characterisation of underlying inflammation.

There is evidence that VOCs can be used to classify patients by asthma phenotype with high accuracy (eg 15-17). VOCs will therefore allow different phenotypes of asthma defined by different inflammatory responses to be distinguished. This clearly points towards the use of breath biopsy as a companion diagnostic for therapy stratification. There is already evidence showing that VOCs can be used to predict steroid responsiveness in mild to moderate asthma sufferers (18). In this study, VOCs were considerably more accurate at predicting responsiveness to steroid treatment (AUC 0.88 ± 0.16) than either F_eNO or measurements of eosinophil cells from sputum samples (see Figure 4).

Personalised Medicine

A reliable tool to tailor treatment to asthma phenotype would be an important development. The past decade has seen the approval of biological drugs targeted at specific inflammatory subtypes, yet their approval for clinical use has been delayed due to their high costs and the difficulty of identifying patients with the correct asthma phenotype who would benefit.

The lack of stratifying diagnostics means that current guidelines advocate a ‘trial and error’ approach, resulting in increased healthcare costs, prolonged periods of poor disease control and an increased risk of exacerbations.

Using biomarkers to characterise asthma subtypes has so far not found widespread adoption due to the invasive nature of existing cell-harvesting techniques, such as bronchoalveolar lavage and the limited applicability of F_eNO testing at identifying responders to new targeted biologic treatments aimed at patients who remain uncontrolled under traditional asthma therapies. The blood-based biomarker IgE

has shown positive predictive potential for treatment response to the targeted anti-IgE treatment XOLAIR®/Omalizumab (19). Therapy failure still occurs in a high percentage of adult and paediatric patients however, and a significant fraction of sufferers lose responsiveness to XOLAIR® over time (20). Additionally, a large portion of patients that would benefit are currently not eligible based on the IgE entry criteria. Similar issues have arisen with the use of periostin and blood eosinophils to identify potential responders to novel asthma therapies.

Promising initial results from the U-BIOPRED consortium project show that VOCs measured by a FAIMS device could be used to stratify asthmatic patients into treatment subgroups. For example, VOCs discriminated between anti-IgE-treated XOLAIR®/Omalizumab and non-treated severe asthma patients with 83% accuracy (21,22).

Future Outlook

VOCs in exhaled breath have enormous potential as biomarkers for disease diagnosis, phenotyping and monitoring. Promising results in the field of asthma indicate that they reflect specific immunological processes that underlie chronic inflammation, making them ideal biomarkers for personalised medicine and treatment stratification.

Highly powered studies such as the MRC EMBER, LuCID and STRATA programmes will generate key evidence for using VOC biomarkers in diagnostic tests (23,8,10). The standardisation of breath sampling technology by the Breathe Free Consortium and the development of analytical techniques like FAIMS with the potential for deployment at point-of-care will allow the widespread adoption of breath biopsy in many clinical applications (4,6).

References

1. van der Schee, MP *et al*, Breathomics in lung disease, *Chest* 147(1): pp224-231, 2015
2. Herbig J and Beauchamp J, Towards standardization in the analysis of breath gas volatiles, *J Breath Res* 8(3): 37101, 2014
3. Risby TH, Critical issues for breath analysis, *J Breath Res* 2: 30302, 2008
4. Visit: www.breathe-free.org
5. Beauchamp J *et al*, On the use of Tedlar® bags for breath-gas sampling and analysis, *J Breath Res* 2(4): 46001, 2008
6. Covington JA *et al*, The application of FAIMS gas analysis in medical diagnostics, *Analyst* 140: pp6,775-6,781, 2015
7. Wilks A *et al*, Characterization of a miniature, ultra-high-field, ion mobility spectrometer, *Int J Ion Mobil Spectrom* 15(3): pp199-222, 2012
8. Visit: www.clinicaltrials.gov/ct2/show/nct02612532
9. Visit: www.owlstonemedical.com/clinical-pipeline/intercept
10. Visit: www.owlstonemedical.com/clinical-pipeline/strata
11. Visit: www.globalasthmareport.org/resources/global_asthma_report_2014.pdf
12. Montuschi P *et al*, Diagnostic performance of an electronic nose, fractional exhaled nitric oxide, and lung function testing in asthma, *Chest* 137(4): pp790-796, 2010
13. Dallinga JW *et al*, Volatile organic compounds in exhaled breath as a diagnostic tool for asthma in children, *Clin Exp Allergy* 40(1): pp68-76, 2010
14. Dragonieri S *et al*, An electronic nose in the discrimination of patients with asthma and controls, *J Allergy Clin Immunol* 120(4): pp856-62, 2007
15. Ibrahim B *et al*, Non-invasive phenotyping using exhaled volatile organic compounds in asthma, *Thorax* 66(9): pp804-809, 2011
16. Fens N *et al*, Exhaled air molecular profiling in relation to inflammatory subtype and activity in COPD, *Eur Respir J* 38(6): pp1,301-1,309, 2011
17. Schleich FN *et al*, Volatile organic compounds discriminate between eosinophilic and neutrophilic inflammation in vitro, *J Breath Res* 10(1): 16006, 2016
18. van der Schee MP *et al*, Predicting steroid responsiveness in patients with asthma using exhaled breath profiling, *Clin Exp Allergy* 43(11): pp1,217-1,225, 2013
19. Westerhof GA *et al*, Biomarkers to identify sputum eosinophilia in different adult asthma phenotypes, *Eur Respir J* 46(3): pp688-696, 2015
20. Visit: www.nice.org.uk/guidance/ta278?unlid=95300737020161031161546
21. Santini G *et al*, Discrimination between oral corticosteroid-treated and oral corticosteroid-non-treated severe asthma patients by an electronic nose platform, *Eur Respir J* 44(58), 2014
22. Santini G *et al*, Breathomics can discriminate between anti IgE-treated and non-treated severe asthma adults, *Eur Respir J* 46, 2015
23. Visit: www.mrc.ac.uk/research/initiatives/stratified-medicine/molecular-pathology/research



Dr Marc van der Schee is Lead Clinical Research Scientist at Owlstone Medical and has nearly 10 years' experience in the field of volatile biomarker research. He holds a medical degree as well as a doctorate in Biomedical Sciences, and obtained a PhD by studying the use of volatile biomarkers for disease diagnosis, monitoring and prognosis prediction. During his PhD, Marc received various awards including a Marie Curie Fellowship.
Email: marc.vanderschee@owlstone.co.uk



Dr Theodore W Wilson is a Scientific Writer at Owlstone Medical. He completed a PhD in Physical Chemistry in 2012, followed by a postdoctoral position in Aerosol Science at the University of Leeds, UK. In 2016, Theodore was awarded a British Science Association Media Fellowship, during which he wrote for the technology section of *The Mirror*.
Email: theo.wilson@owlstone.co.uk

Breath Biopsy® Products

Bring VOC Biomarker research to your own lab with our complete range of products for breath collection and analysis.

Breath collection



ReCIVA™ Breath Sampler

ReCIVA provides an out-of-the-box means to collect breath samples reliably and reproducibly.

- Collect different breath fractions (enriched bronchial or end-tidal)
- Pre-concentrate a large volume of breath onto sorbent tubes, enabling even VOCs present at low levels to be detected
- Developed in conjunction with a multi-disciplinary group of 100 breath analysis experts

owlstonemedical.com/reciva




CASPER Clean Air Supply

CASPER provides a consistent, portable supply of clean air, ideal for use with the ReCIVA™ Breath Sampler.

- Ensure breath samples are not contaminated by external VOCs
- Remove the need for a clean air supply line
- Portable - to wherever the patient is in clinic

owlstonemedical.com/casper

VOC analysis



Lonestar VOC Analyzer

The Lonestar VOC Analyzer is an easy to use instrument for the detection of VOC biomarkers in clinical samples such as breath and other biospecimens.

- Rapid, sensitive, and selective VOC analyzer
- Based on proven FAIMS technology
- Biomarker discovery and targeted analysis

owlstonemedical.com/lonestar

"We have found FAIMS technology able to non-invasively detect a broad range of diseases with high sensitivity and specificity."

Prof. J. Covington, Warwick University





Breath Biopsy[®] Biomarker Services

Discover novel non-invasive biomarkers for cancer, inflammatory and infectious disease with our comprehensive Breath Biomarker Services.

Overview



Owlstone Medical's Biomarker Services enable academics, clinicians and pharmaceutical companies to identify novel breath biomarkers for diagnostics as well as precision medicine applications including therapeutic response, patient stratification and outcome prediction.

We offer a complete set of services for analysis of volatile organic compound (VOC) biomarkers in breath and other biospecimens. Our services include everything from consultation on optimal study design to biomarker discovery and external blinded validation of the test in the context of the clinical use case.

Breath collection

We ensure that breath samples are collected reliably and reproducibly by installing the ReCIVA[™] Breath Sampler at clinical sites, and providing training on optimal breath collection. Breath samples are pre-concentrated on stable sorbent tubes and shipped to Owlstone Medical's service laboratory for analysis.

Biomarker discovery

Owlstone Medical offers the largest laboratory for breath analysis. VOC Biomarkers are analyzed using our proven FAIMS technology. Our analysis experts develop robust classification algorithms from rich spectral FAIMS data. Our exceptional analytical platforms, rigorous quality standards and experienced staff deliver reliable, reproducible and rapid results.



In silico verification and validation

Classification algorithms are verified *in silico* using our proprietary database containing thousands of patient VOC profiles, to provide early evidence of likely clinical performance. Finally, an external blinded test validation is performed in patient groups representative of the intended use population.



Biomarker Services

Identify novel non-invasive
breath biomarkers

owlstonemedical.com/services

