Development of a Compact, IoT-enabled **Electronic Nose for Breath Analysis**

A. Tiele¹, A. Wicaksono¹, S.K. Ayyala¹ and J.A. Covington^{1,*}

¹ School of Engineering, University of Warwick, Coventry, CV4 7AL, UK.

* Corresponding author: Prof. James Covington, <u>J.A.Covington@warwick.ac.uk</u>, +44 (0) 24 7657 4494



Modern breath analysis focuses on the detection of exhaled volatile organic compounds (VOCs) for the diagnosis and/or monitoring of disease. Breath VOCs are the by-products of normal metabolic activity, and in some cases,



ANALYSIS & RESULTS

The device has been tested with chemical standards (acetone, isopropanol and 1-propanol). The PCA results for these tests are shown in Figure 3.

WARWICK THE UNIVERSITY OF WARWICK

specific biomarkers associated with disease (e.g. cancer and COPD [1,2]).

- It has long been suggested that applications of diagnostic breath analysis must extend beyond laboratories and pilot studies to standard clinical practice and home-use [3]. The latter requires a compact, personal and portable diagnostic device, capable of sampling and analysing breath at any time or place.
- Most breath analysis technologies are "high-end" (expensive and complex). However, the latest generation of MEMS technology metal oxide (MOX) gas sensors could provide an alternative low-cost solution.
- While no single analysis technique can provide complete diagnosis of an individual, electronic nose (E-nose) technology has the advantages of being relatively low-cost, low-power, user-friendly and portable.
- In this work, we report on the design and development of a compact, internet-of-things (IoT) enabled E-nose for breath analysis.

MEMS MOX-BASED E-NOSE

The developed unit includes an array of 10 MEMS MOX-based gas sensors (Table 1), including many of the most relevant sensors currently available.



- Exhaled breath samples were collected from 18 volunteers (15 males, 3) females, ages 23-28). A typical output response is shown in Figure 4.
- To simulate a control vs. disease group case-control study, a peppermint breath test was conducted. Subjects provided breath samples pre- and postconsumption of a peppermint oil capsule. This produces a well-defined, but temporary, change in breath composition.
- Classification analysis demonstrates an AUC ± 95% = 0.80, sensitivity = 0.83 (0.59 – 0.96), specificity = 0.72 (0.47 – 0.90), p-value = 0.0007. These results were achieved using the Support Vector Machine (SVM) classifier.

Our system (Figure 1) is compact and uses a microcontroller with Wi-Fi capabilities (ESP32, Espressif Systems, China) for integration with future IoT infrastructure.



Figure 1: Warwick Breath E-nose



Table 1: Deployed Sensors

BREATH SAMPLING

3

Our E-nose design (Figure 2) includes an integrated sampling system, using a heated (35°C) sampling tube (50mL volume) with a $CO_2/Temp/Hum$ gas sensor. This method samples end-tidal breath using the displacement principle.



Figure 5: Radar Plot of Normalised Features

Figure 6: ROC Curve

A radar plot of normalised features is shown in Figure 5 and a generated ROC curve is shown in Figure 6. These results demonstrate that the developed unit is capable of detected changes in exhaled breath.





Figure 2: System Diagram

The functionality of the developed device was demonstrated with the testing of chemical standards and a simplified case-control study using peppermint oil. It is our intention to deploy this system in a UK hospital in upcoming breath research studies.

► REFERENCES

[1] Westhoff, M. et al. Ion mobility spectrometry for the detection of VOCs in exhaled breath of patients with lung cancer: Results of a pilot study. *Thorax*. 2009, 64, 744.

[2] Besa, V. et al. Exhaled volatile organic compounds discriminate patients with chronic obstructive pulmonary disease from healthy subjects. Intl. J. of COPD. 2015, 10, 399-406.

[3] Hunter G.W. and Dweik R.A. Applied breath analysis: an overview of the challenges and opportunities in developing and testing sensor technology for human health monitoring in aerospace and clinical applications. J. Breath Res., vol. 2, no. 3, p. 37020, 2008.