

Determination of N-methylpyrrolidine in cefepime hydrochloride using field asymmetric waveform ion mobility spectrometry combined with mass spectrometry

Colin S Creaser^{1*}, Robert W Smith¹, Lisa Cox¹, James C Reynolds¹, Mark Powell²

Centre for Analytical Science, Department of Chemistry, Loughborough University, Leicestershire, LE11 3TU; ²Quay Pharmaceuticals Ltd., Deeside Industrial Park, Flintshire. CH5 2NS.

c.s.creasr@lboro.ac.uk

Introduction and Background

- Cefepime is a cephalosporin antibiotic with low toxicity and broad range activity.
- N-methyl pyrrolidine (NMP),^{1,2} may be present as an impurity in Cefepime and must not exceed 0.3 %, because of its toxicity.³
- The determination of NMP in Cefepime presents a challenge because, (i) Cefepime rapidly degrades to form NMP in aqueous solution during the analysis, particularly outside the pH range 4-6, and above 4 °C, and (ii) the protonated Cefepime generated by ESI readily fragments in the mass spectrometer interface to form protonated NMP by 'in-source' collision-induced dissociation (in-source CID).
- The hyphenation of high-field asymmetric waveform ion mobility spectrometry with mass spectrometry (FAIMS-MS) allows ions to be separated rapidly on the basis of differential ion mobility and m/z .⁴⁻⁶
- The determination of NMP in Cefepime using FAIMS-MS by the direct infusion of aqueous samples into the ESI source is demonstrated and the quantitative performance evaluated.




Figure 1. Structures of cefepime and impurity NMP

Instrumentation and Experimental

Instrumentation:

- A chip-based prototype FAIMS device (Owlstone Ltd.) was inserted in front of the capillary inlet of an orthogonal acceleration TOF MS (Agilent Technologies 6230 TOF) fitted with a Jet Stream ESI source (Figure 2).
- The FAIMS device has parallel electrode pairs with a 100 μm electrode gap and 700 μm path length. Devices with trench lengths of 78.1 mm and 97.0 mm were; the latter for the quantitative studies.




Figure 2. (a) Photograph of the miniaturised chip-based FAIMS device and (b) Schematic diagram of FAIMS interfaced with the TOF-MS.

Experimental:

- Solutions of NMP and Cefepime prepared in methanol/water (50:50) + 0.1% formic acid were directly infused (10 $\mu\text{L}/\text{min}$) into the ESI source, which was operated in positive ion mode
- ESI-FAIMS-MS conditions: sheath and drying gas flows, 9 and 7 L/min; sheath and drying gas temperatures, 250°C and 150°C; nebuliser pressure, 35 psig; spray shield, 2500 V; transfer capillary voltage and nozzle voltage, 3000 V and 2000 V respectively
- Dispersion fields (DF) of 180 - 300 Td were used for optimising separation and DF 180 Td used for quantitative studies.
- The 'fragmentor' voltage applied in the TOF inlet was varied between 75 and 375 V (25 V steps), but set to 200 V for quantification measurements. The MS acquisition scan rate was 10 scans/s for FAIMS scanning and 1 scan/s for quantification with static FAIMS

Results: In-source CID Fragmentation Studies

- The effect of fragmentor voltage on the intensities of the NMP (m/z 86.07) and cefepime (m/z 481.13) ions was investigated, to determine the extent of cefepime fragmentation by in-source CID.
- NMP and cefepime were infused separately and the fragmentor voltage varied in the range 75-375 V (Figure 3).
- The maximum response for transmission of the NMP ion was at a fragmentor voltage of 200 V (Figure 3a). However, infusion of a Cefepime standard also resulted in a response for NMP as a result of in-source CID at all fragmentor voltages in the range 75-375 V.
- The NMP response in the Cefepime sample directly infused into the ESI source arises from a combination of NMP in the Cefepime sample and in-source CID, preventing the direct determination of NMP.




Figure 3. a plot of intensity of m/z 86.07 (blue) and 481.13 (red) in (a) NMP standard; and (b) Cefepime sample without NMP added

Acknowledgements

We thank Owlstone Ltd and Agilent Technologies for financial and technical support and for the provision of instrumentation, and Quay Pharmaceuticals for providing samples.

Results: ESI-FAIMS-MS of NMP and Cefepime

- Figure 4 shows the mass-selected ion responses for Cefepime (m/z 481) and NMP (m/z 86) in the FAIMS compensation field (CF) spectrum of a solution of Cefepime containing 0.1 % (w/w) NMP.
- The NMP selected ion response (blue trace) shows two peaks, at CF \sim 0 Td, arising from NMP in the Cefepime sample, and at CF 0.7 Td, arising from the transmission of intact Cefepime, followed by in-source CID to form NMP.
- Cefepime and NMP were resolved by FAIMS at a DF of 180 Td (Figure 4, insert).
- Singly and doubly protonated cefepime are the main peaks observed in the mass spectrum of a Cefepime solution without FAIMS separation, (Figure 5a)
- The NMP ion is the base peak in the mass spectrum obtained with FAIMS selection (DF 180 Td, CF 0 Td) (Figure 5b).




Figure 4. FAIMS CF spectrum of NMP (m/z 86.07, blue trace) and cefepime (m/z 481.13, red trace) for a Cefepime sample spiked with NMP (0.1% w/w) at a DF of 180 Td. Insert, the spectrum in the CF range -0.5 to +0.5 Td.

- FAIMS selection allows NMP in Cefepime to be separated from NMP formed by in-source CID.




Figure 5. Mass spectra of NMP (0.05 $\mu\text{g}/\text{ml}$) in cefepime (50 $\mu\text{g}/\text{ml}$): (a) without FAIMS separation and (b) with the FAIMS set to CF 0 Td, DF 180 Td to transmit NMP.

Results: Determination of NMP in Cefepime by ESI-FAIMS-MS

- The FAIMS was operated in 'static' mode (CF 0 Td, DF 180 Td), so that NMP was selectively transmitted through the device for the quantification of NMP in Cefepime. The quantitative performance of direct ESI-FAIMS-MS is summarised in Table 1.
- The LOQ was 0.013 $\mu\text{g}/\text{ml}$, equivalent to 0.027% (w/w) in the Cefepime sample, which is well below the 0.3% threshold
- Good linearity and precision were observed for NMP making direct infusion ESI-FAIMS-MS suitable for the determination of NMP in Cefepime

Table 1. Quantitative performance for the ESI-FAIMS-MS determination of NMP (0.1% w/w) in cefepime (n = 6)

LOD / $\mu\text{g}.\text{ml}^{-1}$	LOQ / $\mu\text{g}.\text{ml}^{-1}$	LDR / $\mu\text{g}.\text{ml}^{-1}$	R ²	%RSD
0.004	0.013	0.005-0.5	0.9979	3.85

Conclusions

- Solution degradation and in-source CID of Cefepime to NMP presents an analytical challenge for the determination of NMP in Cefepime.
- The FAIMS separation of NMP and Cefepime allows NMP in Cefepime to be distinguished from NMP ions generated by in-source CID.
- Direct infusion of a cefepime solution using ESI-FAIMS-MS is a rapid, direct ambient ionization method for quantifying NMP in Cefepime.
- Quantitative determination by ESI-FAIMS-MS was found to be sensitive and precise, with an LOQ well below the required 0.3 % threshold.

References

- Sader, H.S.; Jones, R.N.; *Antimicrob. Newsl.* **1992**, 8, 75-84
- Prasana, J.S.; Sharma, H.K.; Mukkanti, K.; Kumar, V.J.; Raja, G.; Sivakumaran, M.J.; *J. Chrom. Sci.*; **2010**, 48, 830-834
- Farajzadeh, M.A.; Goushuni, L.; Bashour, Y.; *J. Sep. Sci.*; **2010**, 33, 3767-3773
- Brown, L.J.; Smith, R.W. et al. *Anal. Chem.*; **2012**, 84, 4095
- Smith, R.W. et al.; *J. Chrom. A*; **2013**, 1278, 76
- Smith R.W.; Reynolds, J.C.; Lee, S.-L.; Creaser, C.S.; *Anal. Methods*; **2013**, 5, 3799